

Air Showers are self contained chambers installed at entrances to cleanrooms and other controlled environments. They minimize particulate matter entering or exiting the clean space.

Personnel and materials entering or exiting the controlled environment are "scrubbed" by high velocity HEPA-filtered air jets with velocities of 18-30 m/s (3,543 -5,906 fpm). Contaminated air is then drawn through the base within the unit, filtered and recirculated.


Esco is the pioneering company for the construction of air showers with different applications on areas in the micro-electronics, semiconductors, pharmaceutical, spray painting, laboratory animal research and food markets.

#### **Main Features**

- High velocity shower jets in excess of 18-30 m/s (3,543 -5,906 fpm) to ensure efficient scrubbing action to remove particulate matter.
- Operating modes can be programmed in the field.
- Microprocessor controller supervises all functions.
- Mini-pleated H14 filtration achieves > 99.9999% typical efficiency at 0.3 micron particles.
- A disposable prefilter with 85% arrestance extends the life of the main filter.
- An emergency stop button is mounted on both sides of the shower.
- Indicator lights mounted on both sides of the exterior of the air shower to regulate traffic flow, in and out of the cleanroom.
- Permanently lubricated direct drive centrifugal blowers used in conjunction with stainless steel air nozzles.

## **EAS-EUS Airflow System**

- Air is forced by the blower(s) through H14 filter(s) which are 99.9999% efficient against particles of 0.3 microns.
- Filtered air is ejected through nozzles at high velocities into the chamber. These turbulent air streams disperse particulate matter on all surfaces.
- Dispersed particulate matter migrate with the air stream towards the lower areas in the air shower chamber. Air enters the blower supply plenum through pre-filter(s) installed at the base of the chamber.
- The air is continuously filtered and recirculated. The air shower is a selfcontained device and does not exchange air with the environment it is placed in.




H14-filtered air

Unfiltered/potentially contaminated air









#### One-Way


Personnel may enter the controlled environment but not exit through the air shower. At the idle state, the clean side door is locked while the grey side is unlocked. This mode of operation is useful for controlling traffic patterns into the controlled environment.

#### Two-Way One-Way


Personnel may enter or exit the controlled environment through the air shower. The air shower program is able to detect if the person is entering or exiting the controlled environment via door sensors and a time-sequenced control.

#### Two-Way


Personnel may enter or exit the controlled environment through the air shower. In both directions the air shower is activated. This mode of operation is useful in pharmaceutical and lab animal research applications to prevent the egress of hazardous substances and allergens from the controlled environment.



## **Optional EAS Exit Directions**



## **EAS Blowing Type**



Double Side Blowing (EAS-B and C series)

Cleanroom Air Shower • EAS-EUS


per client's specifications


2

	Gei	neral Specifica	ations, Cleanroom Air Show	ver, Model EAS-EUS (A-Serie	es)
	No	te to customer: Ins	sert electrical voltage number into last	t model number digit _ when ordering	g.
Model			EAS-1A_	EAS-2A_	EAS-3A_
External Dimensions (W x D x H)			1385 x 1040 x 2205 mm 54.53" x 40.94" x 86.81"	1385 x 2080 x 2205 mm 54.53" x 81.89" x 86.81"	1385 x 3120 x 2205 mm 54.53" x 122.83" x 86.81"
Internal Work Area, Dimensions (W x D x H)			915 x 920 x 2080 mm 36.02" x 36.22" x 81.89"	915 x 1960 x 2080 mm 36.02" x 77.17" x 81.89"	915 x 3000 x 2080 mm 36.02" x 118.11" x 81.89"
Clear Door Openir	ng (WxH)			815 x 2030 mm (32" x 79.9")	
Air Change (30-18 m/s)		110 - 182/ Hr	105 - 175/ Hr	104 - 172/ Hr	
Initial Airflow Velo	ocity			18-30 m/s (3,543 - 5,906 fpm)	
Number of Nozzles			6	12	18
Air Shower Durati	on			Factory set at 12 seconds (adjustable	)
Persons Per Cycle			1 2-3		
			3 6-9		
Personnel Flow (Persons / Min. )		Above figures based on: Total Cycle Time of 20 seconds (12 seconds of Air Shower + 4 seconds for buffer time + 5 second blower stopping time delay / personnel entrance and exit)			
Filtration Efficiency					
			Pre-Filter: Arrestance 85%, efficiency 20%		
			Main Filter: H14 filter		
Filtration Elements	5		Pre-Filter: Disposable and non-washable polyester fibers		
LED light			1 x 19W 2 x 19W 3 x 19W		
Air Shower Construction		1.5 mm/ 0.06"/18 electro-galvanised steel / White oven-baked epoxy antimicrobial powder-coated finish			
Flectrical -	240-220V, AC, 6	60/50Hz, 1Ø	EAS-1A8	EAS-2A8	EAS-3A8
	120-110V, AC, 60/50Hz, 1Ø		EAS-1A9	EAS-2A9	EAS-3A9
	Assembled (W x D x H)	Pallet A	1680 x 1300 x 2400 mm 66.1" x 51.2" x 94.5"	N/A	N/A
		Pallet B	N/A	N/A	N/A
		Pallet C	N/A	N/A	N/A
		Shipment Weight	Pallet A: 529.2 kg (1166.7 lbs)	N/A	N/A
	Modular (W x D x H)	Pallet A	N/A	1680 x 1300 x 2400 mm 66.1" x 51.2" x 94.5"	1680 x 1300 x 2400 mm 66.1" x 51.2" x 94.5"
Shipping Weight		Pallet B	N/A	1680 x 1300 x 2400 mm 66.1" x 51.2" x 94.5"	1680 x 1300 x 2400 mm 66.1" x 51.2" x 94.5"
and Dimensions (W x D x H)		Pallet C	N/A	N/A	1680 x 1300 x 2400 mm 66.1" x 51.2" x 94.5"
		Shipment Weight	N/A	Pallet A: 449.8 kg (991.6 lbs) Pallet B:449.8 kg (991.6 lbs)	Pallet A: 449.8 kg (991.6 lbs) Pallet B: 371 kg (817.9 lbs) Pallet C: 449.8 kg (991.6 lbs)
	Unassembled (W x D x H)	Pallet A	2500 x 1500 x 1090 mm 98.4" x 59.1" x 42.9"	2500 x 1500 x 1390 mm 98.4" x 59.1" x 54.7"	2500 x 1500 x 1090 mm 98.4" x 59.1" x 42.9"
		Pallet B	N/A	N/A	2500 x 1500 x 1235 mm 98.4" x 59.1" x 48.6"
		Shipment Weight	Pallet A: 530 kg (1168.5 lbs)	Pallet A: 842 kg (1856.3 lbs)	Pallet A: 568 kg (1252.2 lbs) Pallet B:645 kg (1421.9 lbs)
	Assembled		Pallet A: 5.24 m³ (185 ft³)	N/A	N/A
Shipping Volume, Maximum	Modular		N/A	Pallet A: 5.24 m³ (185 ft³) Pallet B: 5.24 m³ (185 ft³)	Pallet A: 5.24 m³ (185 ft³) Pallet B: 5.24 m³ (185 ft³) Pallet C: 5.24 m³ (185 ft³)
	Unassembled		Pallet A: 4.09 m³ (144.4 ft³)	Pallet A: 5.21 m³ (183.9 ft³)	Pallet A: 5.24 m³ (185 ft³) Pallet B: 4.63 m³ (163.5 ft³)

## Esco Cleanroom Air Showers, Modes of Shipment, Model EAS-2A\_-EUS











# **ENGINEERING DRAWING MODEL EAS-EUS (A-SERIES)**

EAS-1A\_-EUS




**EAS-2A\_-EUS** 




EAS-3A\_-EUS



4






- 1. LED Housing
- 2. Esco Sentinel™ Silver Microprocessor Controller
- 3. Unit Circuit Breaker Button
- 4. Unit On/Off Switch
- 5. Stainless Steel Nozzles
- 6. G4 Prefilters (Disposable)
- 7. Electrical Panel
- 8. H14 filters
- 9. Blower
- 10. Indicator Lights (for both grey and clean sides)
- 11. Emergency Stop Button (for both grey and clean sides)


	Ger	neral Specifica	ntions, Cleanroom Air Show	er, Model EAS-EUS (B-Serie	s)*	
	No	ote to customer: In:	sert electrical voltage number into last	model number digit _ when ordering	g.	
Model			EAS-1B_	EAS-2B_	EAS-3B_	
External Dimensions (W x D x H)			1655 x 1040 x 2205 mm 65.16" x 40.94" x 86.81"	1655 x 2080 x 2205 mm 65.16" x 81.89" x 86.81"	1655 x 3120 x 2205 mm 65.16" x 122.83" x 86.81"	
Internal Work Area, Dimensions (W x D x H)			915 x 920 x 2080 mm 36.02" x 36.22" x 81.89"	915 x 1960 x 2080 mm 36.02" x 77.17" x 81.89"	915 x 3000 x 2080 mm 36.02" x 118.11" x 81.89"	
Clear Door Openin	ng (WxH)			815 x 2030 mm (32" x 79.9")		
Air Change (30-18 m/s)		219 - 364/ Hr	209 - 349/ Hr	207 - 344/ Hr		
Initial Airflow Velocity		18-30 m/s (3,543 - 5,906 fpm)				
Number of Nozzles			12	24	36	
Air Shower Duration			Factory set at 12 seconds (adjustable)			
Persons Per Cycle			1	2-3		
			3 6-9			
Personnel Flow (Persons / Min. )		Above figures based on: Total Cycle Time of 20 seconds (12 seconds of Air Shower + 4 seconds for buffer time + 5 second blower stopping time delay / personnel entrance and exit)				
Filtration Efficiency			Main Filter: >99.9999% at 0.3 μm			
			Pre-Filter: Arrestance 85%, efficiency 20%			
			Main Filter: H14 filter			
Filtration Elements			Pre-Filter: Disposable and non-washable polyester fibers			
Fluorescent Lamp			1 x 19W	2 x 19W	3 x 19W	
Air Shower Construction		1.5 mm/ 0.06"/18 electro-galvanised steel / White oven-baked epoxy antimicrobial powder-coated finish				
Electrical	240-220V, AC, 60/50Hz, 1Ø		EAS-1B8	EAS-2B8	EAS-3B8	
	120-110V, AC, 60/50Hz, 1Ø		EAS-1B9	EAS-2B9	EAS-3B9	
	Assembled (W x D x H)	Pallet A	1750 x 1250 2400 mm 68.9" x 49.2" x 94.5"	N/A	N/A	
		Pallet B	N/A	N/A	N/A	
		Pallet C	N/A	N/A	N/A	
		Shipment Weight	Pallet A: 530.5 kg (1169.6 lbs)	N/A	N/A	
	Modular (W x D x H)	Pallet A	N/A	1750 x 1250 x 2400 mm 68.9" x 49.2" x 94.5"	1750 x 1250 x 2400 mm 68.9" x 49.2" x 94.5"	
Shipping Weight and Dimensions (W x D x H)		Pallet B	N/A	1750 x 1250 x 2400 mm 68.9" x 49.2" x 94.5"	1750 x 1250 x 2400 mm 68.9" x 49.2" x 94.5"	
		Pallet C	N/A	N/A	1750 x 1250 x 2400 mm 68.9" x 49.2" x 94.5"	
		Shipment Weight	N/A	Pallet A: 530.5 kg (1169.6 lbs) Pallet B: 530.5 kg (1169.6 lbs)	Pallet A: 530.5 kg (1169.6 lbs) Pallet B: 462 kg (1018.5 lbs) Pallet C: 530.5 kg (1169.6 lbs)	
	Unassembled (W x D x H)	Pallet A	2500 x 1500 x 1090 mm 98.4" x 59.1" x 42.9"	2500 x 1500 x 1090 mm 98.4" x 59.1" x 42.9"	2500 x 1500 x 1380 mm 98.4" x 59.1" x 54.3"	
		Pallet B	N/A	2500 x 1500 x 1090 mm 98.4" x 59.1" x 42.9"	2500 x 1500 x 1455 mm 98.4" x 59.1" x 57.3"	
		Shipment Weight	Pallet A: 608 kg (1340.4 lbs)	Pallet A: 608 kg (1340.4 lbs) Pallet B:452 kg (996.5 lbs)	Pallet A: 730.7 kg (1610.9 lbs) Pallet B:736.7 kg (1624.1 lbs)	
	Assembled		Pallet A: 5.25 m³ (185.4 ft³)	N/A	N/A	
Shipping Volume, Maximum	Modular		N/A	Pallet A: 5.25 m³ (185.4 ft³) Pallet B: 5.25 m³ (185.4 ft³)	Pallet A: 5.25 m³ (185.4 ft³) Pallet B: 5.25 m³ (185.4 ft³) Pallet C: 5.25 m³ (185.4 ft³)	
	Unassembled		Pallet A: 4.09 m³ (144.4 ft³)	Pallet A: 4.09 m³ (144.4 ft³) Pallet B: 4.09 m³ (144.4 ft³)	Pallet A: 5.18 m³ (182.9 ft³) Pallet B: 5.46 m³ (192.8 ft³)	

Note: EAS-B series have an option to have NO BASE with guide rail on the side of each module.

## Esco Cleanroom Air Showers, Modes of Shipment, Model EAS-2B\_-EUS











## **ENGINEERING DRAWING MODEL EAS-EUS (B-SERIES)**


EAS-1B\_-EUS




**EAS-2B\_-EUS** 



EAS-3B\_-EUS



6







- 1. LED Housing
- 2. Esco Sentinel™ Silver Microprocessor Controller
- 3. Unit Circuit Breaker Button
- 4. Unit On/Off Switch
- 5. Stainless Steel Nozzles
- 6. G4 Prefilters (Disposable)
- 7. Electrical Panel
- 8. H14 filters
- 9. Blower
- 10. Indicator Lights (for both grey and clean sides)
- 11. Emergency Stop Button (for both grey and clean sides)

		•	tions, Cleanroom Air Shower, Model EAS			
		e to customer: Ins	ert electrical voltage number into last model number dig			
Model			EAS-2C_	EAS-3C_		
External Dimensions (W x D x H)			2435 x 2080 x 2205 mm 95.87" x 81.89" x 86.81"	2435 x 3120 x 2205 mm 95.87" x 122.83" x 86.81"		
Internal Work Area, Dimensions (W x D x H)			1690 X 1960 X 2080 mm 66.54" x 77.17" x 81.89"	1690 X 3000 X 2080 mm 66.54" x 118.11" x 81.89"		
Clear Door Openin	ıg (WxH)		1590 x 2030 m	ım (62.6" x 79.9")		
Air Change (30-18	m/s)		113 - 188/ Hr	112 - 186/ Hr		
Initial Airflow Velc				543 - 5,906 fpm)		
Number of Nozzle			24	36		
Air Shower Duration	on		Factory set at 12 seconds (adjustable)			
Persons Per Cycle			2-3	4-6		
Personnel Flow (Persons / Min. )			6-9 10-21  Above figures based on: Total Cycle Time of 20 seconds (12 seconds of Air Shower + 4 seconds for buffer time + 5 second blower stopping time delay / personnel entrance and exit)			
Filtration Efficiency						
			Pre-Filter: Arrestance	Pre-Filter: Arrestance 85%, efficiency 20%		
			Main Filter: H14 filter			
Filtration Elements			Pre-Filter: Disposable and r	non-washable polyester fibers		
Fluorescent Lamp			4 x 19W 6 x 19W			
Air Shower Construction			1.5 mm/ 0.06"/18 electro-galvanised steel / White oven-baked epoxy antimicrobial powder-coated finis			
Electrical	240-220V, AC, 60/50Hz, 1Ø		EAS-2C8	EAS-3C8		
	120-110V, AC, 60/50Hz, 1Ø		EAS-2C9	EAS-3C9		
		Pallet A	N/A	N/A		
	Assembled	Pallet B	N/A	N/A		
	(W x D x H)	Pallet C	N/A	N/A		
		Shipment Weight	N/A	N/A		
	Modular	Pallet A	2600 X 1200 X 2400 mm 102.4" x 47.2" x 94.5"	2600 x 1200 x 2400 mm 102.4" x 47.2" x 94.5"		
Shipping Weight		Pallet B	2600 X 1200 X 2400 mm 102.4" x 47.2" x 94.5"	2600 x 1200 x 2400 mm 102.4" x 47.2" x 94.5"		
and Dimensions (W x D x H)	(W x D x H)	Pallet C	N/A	2600 x 1200 x 2400 mm 102.4" x 47.2" x 94.5"		
		Shipment Weight	Pallet A: 753 kg (1660.08 lbs) Pallet B: 753 kg (1660.08 lbs)	Pallet A: 753 kg (1660.08 lbs) Pallet B: 642 kg (1415.4 lbs) Pallet C: 753 kg (1660.08 lbs)		
	Unassembled (W x D x H)	Pallet A	2700 x 1570 x 1890 mm 106.3" x 61.8" x 74.4"	2700 x 1570 x 1890 mm 106.3" x 61.8" x 74.4"		
		Pallet B	2700 x 1570 x 1190 mm 106.3" x 61.8" x 46.9"	2700 x 1570 x 1530 mm 106.3" x 61.8" x 60.2"		
		Shipment Weight	Pallet A: 684.8 kg (1509.7 lbs) Pallet B: 741 kg (1633.63 lbs)	Pallet A: 843 kg (1858.5 lbs) Pallet B: 1073.3 kg (2366.2 lbs)		
	Assembled		N/A N/A			
Shipping Volume, Maximum	Modular		Pallet A: 7.49 m³ (264.5 ft³) Pallet B: 7.49 m³ (264.5 ft³)	Pallet A: 7.49 m³ (264.5 ft³) Pallet B: 7.49 m³ (264.5 ft³) Pallet C: 7.49 m³ (264.5 ft³)		
	Unassembled		Pallet A: 8.01 m³ (282.9 ft³) Pallet B: 5.04 m³ (177.9 ft³)	Pallet A: 8.01 m³ (282.9 ft³) Pallet B: 6.49 m³ (229.2 ft³)		

Note: EAS-C series have an option to have NO BASE with guide rail on the side of each module.

# Esco Cleanroom Air Showers, Modes of Shipment, Model EAS-2C\_-EUS









# **ENGINEERING DRAWING MODEL EAS-EUS (C-SERIES)**

## EAS-2C\_-EUS





## EAS-3C\_-EUS



1. LED Housing

8

- Esco Sentinel™ Silver Microprocessor Controller
- Unit Circuit Breaker Button
- Unit On/Off Switch

- Stainless Steel Nozzles
- G4 Prefilters (Disposable)
- Electrical Panel
- 8. H14 filters
- Blower

- 10. Indicator Lights (for both grey and clean sides)
- 11. Emergency Stop Button (for both grey and clean sides)

#### **Esco Pharma**

21 Changi South Street 1 Singapore 486777 • Tel: +65 65420833 Fax: +65 65426920 • Email: csis.pharma@escoglobal.com

#### Esco Technologies, Inc.

2512 Metropolitan Drive, Suite 120 B Feasterville- Trevose, PA 19053-6738 Tel: 215 322 2155 • Email: eti.pharma@escoglobal.com

#### Esco Gb Ltd

Unit 2 R-evolution @ Gateway 36, Kestrel Way, Barnsley, S70 5SZ Tel: +44 (0) 1226 360 799 • Email: egb.pharma@escoglobal.com



















Esco Global Offices: Bangladesh | Cameroon | China | Denmark | Germany | Hong Kong | Indonesia | Italy | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam